


ADVISORY NOTICE

Products Affected: ExSept Plus and Alcavis 50


Date Issued: 04/16/2007

ADN No.: 001

- ExSept Plus and Alcavis 50 solutions when stored appropriately (see SDS for storage instructions) are stable for 30 months from date of manufacture.
- When opened, ExSept Plus and Alcavis 50 solutions will maintain stability for 180 days, if appropriately stored.
- Frozen bottles if allowed to thaw at room temperature will maintain the labeled concentration. DO NOT warm/heat the bottle to expedite the thawing process.



Quality / Regulatory Director



Date