

Dermacoll-Thin™

WOUND DRESSING

Dermacoll-Thin™ is a gelatin-free hydrocolloid adhesive wound dressing laminated to a semi-permeable, and waterproof polyurethane film backing. Individually packed and sterilized by radiation.

INDICATIONS

- Light to Moderately Exudating Wounds
- Chronic Wounds
- Diabetic Ulcers
- Leg Ulcers
- Pressure Ulcers
- Abrasions



CONTRAINDICATED

For use on third degree burns or exposed muscle, bone and tendons, or infected wounds.

Non adhering to wounds. Dressings can be easily removed without causing pain and without damaging the granulation and epithelial tissue.



ANGELINI

CONTACT US

Angelini Pharma Inc.

8322 Helgerman Court - Gaithersburg - MD 20877 -USA

1 (800) 726-2308


www.Angelini-us.com

Dermacoll-Thin™

WOUND DRESSING

Dermacoll-Thin™ forms a gel upon absorbing wound exudate providing a moist wound environment supporting granulation and epithelialization.

Gel formation prevents Dermacoll-Thin™ from adhering to wounds and the dressing can be removed without causing pain and without damaging the granulation and epithelial tissue.



- Gelatin Free, minimizes associated allergies and odors
- Semi-permeable, waterproof polyurethane backing
- High absorption capacity of Dermacoll-Thin™ binds away excess exudates extending dressing change intervals
- Accelerates the healing process by providing an ideal environment for the formation of granulation tissue
- Slightly transparent, Dermacoll-Thin™ enables wound inspection

Available in the following sizes:

Dermacoll-Thin™ - 4x4 10/Box - **Product Code:** 317044T